

Continuous Run Time 5 hours	☒
5 Adjustable Brightness Degree	☒
5 Color Mode	☒
Inverse Mode	☒
Enhance/Fine Mode	☒
Stand-by Mode	☒
Depth Mode	☒
Vein Image Storage and export	☒
USB port for data transfer	☒
Auto wake up function	☒
Distance sensor and hint	☒
Vein image centre line	☒
Accessory: Enhancer	☒
2 Optional stands: Fixed Stand	☒
Optional stand: Mobile Stand	☒
The device is chargeable in the stand	☒



BLZ Technology (Wuhan) Co., Ltd.

TEL +86-27-6552 0255

WEB www.veinsight.com

Email: info@veinsight.com

ADD: Room 2166, 2176 Building E, Cyberport Dongxing Road, East Lake
New Technology Development Zone, 430074, Wuhan, China

VEINSIGHT

VS500

Lighter, Smaller
and Smarter
Vein Finder



reddot winner 2020



reddot winner 2020
innovative product



BLZ Technology (Wuhan) Co., Ltd.

VEINSIGHT™ vein finder can digitally displays a map of the veins on the surface of the skin in real time, aiding healthcare professionals to verify vein patency. Vein finder is a proven modality that can improve patient care and save valuable time for both clinicians and patients. As one of the most advanced vein finder ever, VS500 has many unique features that original designed by BLZ Tech team to improve the using expirience among end user.


Distance Sensing Technology
Wave Hands to Wake Up VS500

VS500 has distance sensing feature, which will help the end user find perfet distance. You can also wave your hand beneath the device then the device will be automatically activated.




Saving Real Time Vein Image
Exporting to Computer

VS500 can save the vein image in Real time, and conect to an computer to export all saved images.



Venipuncture Point Hint
Provide More Information

VS500 has the unique function called Venipuncture Point Hint mode, which can provide the medical staff the center line of the veins.



Stand is essential in the daily use of vein finder, it can help free both hand for the medical staff.

Application

